

Fixing 3D Model for Printing

CHAT-0180 / REV EXPORT-03 / 2026.02.09

A conversation-specific public summary of fixing 3d model for printing, distilled from the supplied export. The detailed excerpts remain available in the released note file.

STATUS	RECOVERED / SPECIFIC NOTE
CLASSIFICATION	SYS
MATERIAL	CHAT EXPORT / BUILD & TOOLS
DISPOSITION	PUBLIC CONVERSATION EXCERPTS / REDACTED
CATEGORY	Build & Tools
RECORD TYPE	build-thread

TECHNICAL NOTES

Conversation-specific knowledge extracted from the fixing 3d model for printing discussion.

CONVERSATION / QUESTION AND CONTEXT

how do i fi this so i can print it I am using 3ds ma and orca as the slicer

CONVERSATION / EXTRACTED GUIDANCE

-
- An open edge (hole) inside the mouth
- Intersecting surfaces (teeth, lips, inner mouth not merged)
- Multiple shells that are touching but not welded
- Normals flipped in the mouth cavity
- Select the statue
- Right-click -> Convert to Editable Poly
- Go to Modifier List -> ProBoolean

SPECIFICATION

SOURCE	Supplied ChatGPT export
CATEGORY	Build & Tools

CONVERSATION 3 user turns / 4 assistant turns

ATTACHMENTS 14 references detected

PUBLIC DETAIL Extracted knowledge + excerpts

DEVELOPMENT HISTORY

2026.02.09 / EXPORT-04

Matched to the supplied export; key guidance is surfaced separately and the raw transcript remains unpublished.

RELEASED MATERIAL / SOURCE LINKS

CHAT-0180-public-note.txt

PUBLIC-SAFE CONVERSATION NOTE / UNIQUE FILE

<https://u225x.com/downloads/chat-records/CHAT-0180-public-note.txt>

RETURN TO THE LIVING RECORD

Fixing 3D Model for Printing

This PDF preserves the public notes at this revision. Visit the record for updates and attached images.

<https://u225x.com/record/CHAT-0180>